

bioabsorbable atelocollagen hemostatic felt

Cat. No.:	Size
008	10 x 30 mm, 5 pieces

LOT

Batch number

REF

Catalog number



Production date



Expiry date



For single use only

STERILE R

Sterilized by gamma-radiation

CE₁₀₂₃

Notified Body
ITC Zlín, a.s., Czech Republic

GB

Composition: Bovine atelocollagen type I, 99.9% crystalline

Indication group: Haemostatic

Description: Atelocollagen is a very effective haemostatic. In addition to its excellent properties, Hypro-Sorb NT is well tolerated in the tissue. It is non-immunogenic, resorbable and supports the healing processes.

Activity: Hypro-Sorb NT has a specific activity to thrombocytes and releases clotting factors, which together with plasma binding factors promote fibrinogenesis.

Physiological Effects: Hypro-Sorb NT implanted into the tissue resorbs spontaneously in two to four weeks, depending on the amount and the body region.

Indication: Hypro-Sorb NT is indicated for the control of capillary and parenchymatous bleeding. Hypro-Sorb NT is used to stop nasal haemorrhage.

Contraindication: Hypro-Sorb NT is not destined for use in instances of pumping arterial haemorrhage. Hypro-Sorb NT must not be left inside the bone fractures fixed by acrylate adhesives because it considerably reduces the bonding strength.

Side Effects: Hypro-Sorb NT is not recommended for persons known to be hypersensitive to the products of bovine origin. Hypro-Sorb NT should not be left between the skin edges as this could interfere with the wound healing and develop wide scars.

In the case of coagulatory disorders or anticoagulation treatment, the haemostatic effect of Hypro-Sorb NT may be reduced or eliminated.

GB

Warning: Hypro-Sorb NT is sterile. Check original packaging before opening. Do not use if packaging has been opened. Unused Hypro-Sorb NT cannot be re-sterilized. Hypro-Sorb NT is not intended to replace ligation or direct compression in cases of venous bleeding.

Dosage: Hypro-Sorb NT should be used in amounts which fully cover the wounded surface. It is recommended that the surgeon remove the excess felt before closing the wound or leave only the minimum amount in the tissue because it slightly swells and it could exert pressure on the vicinity of the wound. Hypro-Sorb NT implanted into the tissue resorbs spontaneously in two to four weeks depending on its amount and the body region.

Handling instructions: Hypro-Sorb NT is applied usually with its smooth side onto the wound surface, softly pressed down and left on the wound until the fibrine adhesion develops. Bleeding usually stops in two to five minutes; it lasts longer in haemophiliacs. Usually it is used in dry state but it may also be moistened with sterile saline.

Storage conditions: Hypro-Sorb NT must be stored in a dry place and at a normal temperature. Hypro-Sorb NT is not damaged by temperatures from -25°C/-13°F to +50°C/+122°F. It needs to avoid direct sun light.

Producer: Hypro Otrokovice, s. r. o., Přístavní 568, 765 02 Otrokovice, Czech Republic, phone +420 577159727, fax +420 577159724, www.hypro.cz, e-mail: hypro@hypro.cz

Registration: EC Certificate No. 09 0627 QS/NB, EC Design-Examination Certificate No. 09 0628 CN/NB and Declaration of Conformity.

